

Concentricity and Neck Thickness Gauge (1200 R)



Location **Gauteng, Midrand**
<https://www.freeadsz.co.za/x-226993-z>











Measures:

- Fired case run-out (chamber run-out),
- Sized case run-out (sizing dies and press), and
- Loaded round run-out (seating dies and press).
- As well as neck thickness of cases.

Accuracy to 0.01 mm.

Collection in Midrand, or add R99 for Postnet.

R1000 for concentricity gauge only without neck.

 <p>Concentricity and Neck Thickness Gauge</p> <p>https://www.freeadsz.co.za/x-226993-z</p>	 <p>Concentricity and Neck Thickness Gauge</p> <p>https://www.freeadsz.co.za/x-226993-z</p>	 <p>Concentricity and Neck Thickness Gauge</p> <p>https://www.freeadsz.co.za/x-226993-z</p>	 <p>Concentricity and Neck Thickness Gauge</p> <p>https://www.freeadsz.co.za/x-226993-z</p>	 <p>Concentricity and Neck Thickness Gauge</p> <p>https://www.freeadsz.co.za/x-226993-z</p>	 <p>Concentricity and Neck Thickness Gauge</p> <p>https://www.freeadsz.co.za/x-226993-z</p>	 <p>Concentricity and Neck Thickness Gauge</p> <p>https://www.freeadsz.co.za/x-226993-z</p>	 <p>Concentricity and Neck Thickness Gauge</p> <p>https://www.freeadsz.co.za/x-226993-z</p>	 <p>Concentricity and Neck Thickness Gauge</p> <p>https://www.freeadsz.co.za/x-226993-z</p>	 <p>Concentricity and Neck Thickness Gauge</p> <p>https://www.freeadsz.co.za/x-226993-z</p>
---	--	--	--	---	--	--	--	--	--